

Durable nitrile gloves offering excellent comfort and tactile sensitivity across numerous healthcare applications

Improved durability: Owing to Ansell's unique nitrile formulation, MICROFLEX[®] 92-245 exam gloves ensure greater durability, while delivering a comfortable fit and feel

Enhanced dexterity: Ideal for a broad range of healthcare applications, these powder-free nitrile gloves are designed to provide the tactile sensitivity and dexterity needed for finer hand movements

Heightened grip: In addition, MICROFLEX[®] 92-245 medical examination gloves' textured fingertips afford a strong grip, resulting in the secure handling of instruments



Industries

- Healthcare
- Dental
- Veterinary

Recommended For

- General patient care
- General dentistry
- Positioning patients
- Standard patient examinations

TECHNICAL DATA SHEET

| Product Information | |
|------------------------|--|
| Material | Nitrile |
| External Glove Surface | Textured Fingertips |
| Audit Standards | ISO 9001, Medical Device Regulation (EU) 2017/745 Class I |
| Product Certification | Medical Device Regulation (EU) 2017/745 Class I, ASTM D6319, CE Number, EN 455 1-4, IS 15354, ISO 11193, ISO 16604, ISO 9001 |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |
| Country Of Origin | PRC |
| Product Segmentation | Multi-Purpose |
| Antistatic | No |
| Silicone Free | No |

PHYSICAL PROPERTIES

| | Typical Values | Testing Method | |
|---|---------------------|--------------------|------------------|
| Length (mm/inches) | 240 / 9.45 | ASTM D6319 | |
| Freedom from Holes (Inspection level I) | 1.5 AQL | ASTM D5151 | |
| Palm Thickness (mm/mils) | 0.09 / 3.5 | ASTM D6319 | |
| Finger Thickness (mm/mils) | 0.14 / 5.5 | ASTM D6319 | |
| | BEFORE AGING | AFTER AGING | |
| Ultimate Tensile Strength (MPa) | Min 14 | Min 14 | ASTM D412 & D573 |
| Elongation at Break (%) | Min 500 | Min 400 | ASTM D412 & D573 |
| Force at break (N) | Min 7 | Min 6 | EN 455-2 |

ORDERING INFORMATION

| Size | XS (5.5 - 6) | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|--------------|--------------|-------------|-------------|-------------|---------------|
| Product Code | 92245060 | 92134070 | 92134080 | 92134090 | 92245100 |

Performance Standards and Regulatory Compliance



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

MY +603 8310 6688
apac.medical@ansell.com
CN: +86 21 38275000
infochina@ansell.com
JP: +813 5549 8151
info.medical.jp@ansell.com

North America Region

Ansell Healthcare Products LLC
T: +1 800 800 0444
F: +1 800 800 0445

Latin America & Caribbean Region

BR +55 (11) 3356-3100
latam.medical@ansell.com
MX: +52 442 296 2050
latam.medical@ansell.com
CO: +57 1 288 3247
latam.medical@ansell.com

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

Russia

Ansell PYC
Ten. +7 495 258 13 16

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

