



Latex Examination Gloves Powdered , Standard Cuff

MEDIGEN® powdered latex gloves are designed for low to moderate-risk environments, offering effective perspiration absorption that ensures a comfortable and dry feel—even after extended wear.



MEDIGEN® Latex Examination

Length (mm)		
	≥ 230	
Palm (centre of Palm)	0.11 ± 0.03	
Finger (13mm ± 3mm from tip)	0.13 ± 0.03	
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 650	≥ 500
Watertightness	G1	1.5
Physical Dimensions	S2	4.0
Tensile Strength	S2	4.0
Visual Inspection (Major)	S4	2.5
Visual Inspection (Minor)	S4	4.0
Powder Amount	N = 2	≤ 150 mg/glove
Protein Content	N = 3	≤ 200 µg/g

*Full test reports on Sensitisation Testing, Primary Skin Irritation and Residual Proteins Testing are available upon request.

GLOVE SIZES

- S - SMALL
- M - MEDIUM
- L - LARGE

FEATURES

- Smooth finger tips
- Lightly powdered
- Made with natural rubber latex
- Protein level ≤200µg/g
- Ambidextrous
- Natural colour

PACKAGING

- 100 gloves per box
- 10 boxes per carton

REGULATORY COMPLIANCE

- FDA 510(k), MDD 93/42/EEC,
- ARTG 127794, CMDR SOR-98-282,
- REACH, EU 2016/425

STANDARDS

- ISO 11193-1, AS/NZS 4011-1, ASTM D5151,
- ASTM D6124, ASTM D3578, ASTM D5712,
- EN ISO 374-1 (Type B), EN 16523-1,
- EN 374 part 2, 4 & 5, EN 455 part 1, 2 & 3,
- ISO 10993 part 5 & 10, HACCP Certified

MANUFACTURING ACCREDITATIONS

- ISO 9001, ISO 14001, EN ISO 13485

A brand by



For additional information visit us at www.atomsurgical.com or call us at

Atom Surgical
T: +91 9429488976
E: info@atomsurgical.com

Atom Surgical® and related marks are trademarks owned by Atom Surgical or its affiliates. All Rights Reserved.