

DEVICE SPECS

# SURGIFLO™ Haemostatic Matrix



**SURGIFLO™ Haemostatic Matrix with Thrombin**



**SURGIFLO™ Haemostatic Matrix**



**SURGIFLO™ Endoscopic Applicator**

<b>Codes</b>	MS0012	MS0010	MS1995
<b>Description</b>	Intended for haemostatic use by applying to a bleeding surface	Intended for haemostatic use by applying to a bleeding surface	A minimally reflective stainless steel cannula and a stylet (obturator)
<b>Indication</b>	SURGIFLO™ Haemostatic Matrix Kit with Thrombin (SURGIFLO™) is intended for haemostatic use by applying to a bleeding surface	SURGIFLO™ is indicated for surgical procedures (except ophthalmic) for haemostasis, when control of capillary, venous and arteriolar bleeding by pressure, ligature, and other conventional procedures is ineffective or impractical.	The SURGIFLO™ Endoscopic Applicator is intended for use in delivering hemostatic agents to bleeding surgical sites through a 5 mm trocar or larger.
<b>Size</b>	8 mL	8 mL	34cm
<b>Hemostatic Agent &amp; Accessory</b>	Haemostatic Agent	Haemostatic Agent	Accessory
<b>Storage Requirements</b>	Stored dry at controlled room temperature (2°C-25°C)	Stored dry at controlled room temperature (2°C-25°C)	Store at controlled room temperature (15°-30°C)
<b>Preparation Time</b>	Prepare the Thrombin Solution inside the sterile field	Prepare the Matrix inside the sterile field (6 passes to mix)	Ready out of package
<b>Material/Composition</b>	Porcine Gelatin and Human Thrombin	Porcine Gelatin	Stainless Steel Cannula and Plastic Obturator
<b>Qty./Box</b>	6	6	6